



TECHNICAL SPECIFICATIONS

---

**ORTHO VISION® MAX SWIFT ANALYZER**

FOR ID-MTS™ GEL CARDS

# SPECIALLY DESIGNED FOR YOUR LAB. BECAUSE EVERY TEST IS A LIFE.™

## INTENDED USE

ORTHO VISION® Max *Swift* Analyzer automates in vitro Immunohematology testing of human blood utilizing ID-MTS™ Gel Card technology.

### ORTHO VISION Max *Swift* Analyzer:

- Automates test processing functions including liquid pipetting, reagent handling, incubation, centrifugation, reaction grading and interpretation, and data management requirements using ID-MTS Gel Cards and digital image processing.
- Standalone instrument or interfaced to the Laboratory Information System (LIS).
- Fully compatible with ORTHO CONNECT™ to maximize the efficiency of your data, workflow, and resources for integration across your lab or network.

## INSTALLATION AND SITE SPECIFICATIONS

Although trained service personnel install the ORTHO VISION Max *Swift* Analyzer at the laboratory site, the site must be prepared according to site specifications.

## ANALYZER SPECIFICATIONS

### Physical Dimensions

#### ■ ORTHO VISION Max *Swift* Analyzer

Width: 169 cm (66.5 in.)  
Depth: 77 cm (30.3 in.)  
Height: 93.9 cm (37 in.)  
Height with maintenance door open: 133.8 cm (52.7 in.)

#### ■ ORTHO VISION Max *Swift* Analyzer Table Dimensions (highly recommended)

Table without shelf:  
182.8 cm (72 in.) x 76.2 cm (30 in.) x 76.2 cm (30 in.)

Table with shelf at the front:  
182.8 cm (72 in.) x 106.7 cm (42 in.) x 76.2 cm (30 in.)  
Occupies 61 cm (24 in.) in depth on left front of table.

Table with the shelf on the side:  
243.8 cm (96 in.) x 76.2 cm (30 in.) x 76.2 cm (30 in.)  
Occupies 61 cm (24 in.) in width on left side of table.



MAX ANALYZER  
ORTHO VISION®  
*Swift*

1. Monitor
2. Dual Load Station for Samples and Reagents
3. Dual Purpose Drawer
4. ID-MTS Gel Card Waste Drawer
5. Supply Drawer
6. Liquid System
7. Shelf Attachment Positions

## Electrical Data

- Data voltage: 100–240 V AC
- Transient overvoltage: Category II
- Frequency: 50/60 Hz
- Power consumption: 1,000 VA
- Power requirement: One dedicated (3-wired single phase; line-neutral and single circuit) AC power line for connection to facility power.

## Hardware

16 GB RAM, 1 TB Solid State Hard Drive

## Software

- Operating system: Windows 10 IoT Enterprise
- Antivirus system: Cylance Artificial Intelligence Platform (in development)

## Weights

ORTHO VISION Max *Swift* Analyzer: 330 Kg (727.5 lbs)  
ORTHO VISION Max *Swift* Analyzer Table: 113.4 kg (250 lbs)

## Environmental Specifications

- Operating temperature: 18–30° C (64.4–86.0° F)
- Site relative humidity: 15–85% RH (non-condensing)
- Maximum altitude: 2438.0 m (8000 ft.)
- Heat output: 3412 BTU/hr
- Noise level: Average 52 dB(A)

## System Performance

- **Samples:** Up to 84 samples can be loaded, 12 Sample Racks holding 7 tubes each
- **ID-MTS Gel Card Supply Drawer:** 200 ID-MTS Gel Cards (10 sleeves of 20)
- **Reagent Red Blood Cell Supply:** Supports 6 independent RBC racks, each can hold 11-3 mL vials or 6-10 mL vials, including various combinations. The 10 mL rack can accommodate reverse grouping and antibody detection (screen) RRBC. The 3 mL rack can accommodate reverse group, antibody detection and antibody identification RRBC.
- **Diluent:** 4 100 mL positions to accommodate MTS™ Diluent 2 and MTS™ Diluent 2 Plus; 4 10 mL positions to accommodate Titer diluent; Onboard liquid reagent capacity of 4.7 L of saline and 1 L of deionized water; Extended walkaway kit with additional 5.2 L saline capacity available
- **Probes:** 4 pipetting probes to maximize efficiency
- **Heated Incubator:** 12 ID-MTS Gel Cards
- **Room Temperature Incubator:** 24 ID-MTS Gel Cards

- **Centrifuge:** 2 centrifuges; 10 ID-MTS Gel Cards per centrifuge
- **Waste:** 140 ID-MTS Gel Cards; External 10 L liquid waste drain kit available; Optional off-board plumbing draining available

## SAMPLE SPECIFICATIONS

### Supported Sample Types

- Centrifuged whole blood
- Plasma and Serum
- Packed red blood cells

### Supported Sample Tube Sizes

- 16 x 100 mm, 16 x 75 mm
- 2-13 x 100 mm, 12-13 x 75 mm
- 10.25 x 75 mm, 10.25 x 64 mm, 10.25 x 47 mm
- Wide variety of Sarstedt Tubes
- 2.0 mL and 1.5 mL micro-collection containers
- Wide variety of Pediatric Tubes

For additional sample tubes supported please refer to the sample rack selection guide.

### Sample and Test Processing

Continuous, random, STAT access, and batch

### Sample Handling

- Automatic cap detection without workflow interruption
- Automatic clot detection
- Automatic liquid level detection

### Sample and Test Performance Characteristics

Column Agglutination Technology using ID-MTS Gel Cards

#### Testing

ABO/Rh	Direct Antiglobulin
Monoclonal Rh Phenotyping	Antigen Testing
Antibody Screen	QC Testing
Antibody Identification	Selected Cell Panels
Crossmatch	
Serial Dilutions for Titration Studies	

#### Reagent Types

ID-MTS Gel Cards	Diluents
Reagent Red Blood Cells	Quality Control

#### Specimen Types

Plasma
Packed Red Blood Cells
Serum



## CONNECTIVITY

### e-Connectivity® Remote Connectivity

Provides the ability to connect your system to Ortho Clinical Diagnostics in a way that enables remote diagnostics, as well as monitor and review system configuration, data, and performance.

e-Connectivity Technology Ready Network Connection Requirements:

- Category 5e cable with male RJ45 connector for SSL connection included.
- Continuous broadband connection or direct connection to the site LAN with access to the Internet at a speed greater than or equal to 128 kbps.
- Support local area network port speeds of automatic, 10/100 Mbps with half and full duplex.
- Dynamic or Static IP Address, Subnet Mask, and Default Gateway IP Address assigned by the IT department and provided to Ortho.
- Female RJ45 connector on network port within 20 feet of system center.

### Broadband Internet

A broadband internet connection is required.

### System Computer and Interface Specifications

- **Interface Specification:** Bidirectional protocols for a Laboratory Information System (LIS).
- **Remote Review Capable:** An external computer on the laboratory's network where Authorized personnel on an external computer can review results.

### LIS Specifications

- **LIS Interface:** Interface via one of three user configurable physical interfaces.
  - ASTM over RS-232
  - ASTM over TCP/IP
  - Network shared folders
- **ASTM Protocol:** Configurable to one of three options.
  - Basic ASTM (no 'M' record)
  - Enhanced ASTM
  - Vision ASTM, which adds in addition to Enhanced ASTM: Error upload message if the order could not be processed; The LIS can download multi-tube orders; The LIS can send QC orders that specify the card reagent lots to use as well as the expected results
- **Barcode Symbolologies**

NW7 (Codabar)	Code 3 of 9 (Code 39)
ISBT 128	Code 2 of 5 (Interleaved)
Code 128 (A, B, & C subtypes)	

### Communication Ports

- I/O ports include:
  - 1 DB-9 serial port (RS-232 port for LIS support)
  - 1 RJ45 LAN port supports port speeds of automatic, 10/100 Mbps with half and full duplex.
- 5 V 2.0/V 1.1 USB ports are available for printer, handheld barcode reader, other devices

### Printer Specifications

ORTHO VISION Max *Swift* Analyzer can be connected either to a network printer or to a local printer.

Because **Every Test is a Life**™

All trademarks are the property of Ortho Clinical Diagnostics Inc.

Ortho Clinical Diagnostics