

## **Brennan & Co - Products available on HSE 4519**



# Table of Contents

## Products available on HSE 4519

Product	Page
<b>Contacts</b>	<b>3</b>
<b>Actico®</b> - Cohesive Inelastic Bandage	<b>4</b>
<b>Actico® 2C</b> - Two Component Leg Ulcer Kit	<b>4</b>
<b>Actifast®</b> - 2 Way Stretch Tubular Retention	<b>5</b>
<b>Activa®</b> - Hosiery Liners	<b>6</b>
<b>Activa®</b> - Leg Ulcer Hosiery Kits	<b>6</b>
<b>Flexiban®</b> - Padding Bandage	<b>7</b>
<b>Debrisoft®</b> - For Debridement	<b>7</b>
<b>Lomatuell® Pro</b> - Gell Forming Wound Contact Layer	<b>8</b>
<b>ReadyWrap®</b> - The Easy Solution For Self-Care	<b>8</b>
<b>Vliwasorb Adhesive®</b> - Superabsorbent Wound Dressing Compression	<b>9</b>
<b>Vliwasorb® Pro</b> - Highly Comfortable Dressing For Highly Exuding Wounds	<b>9</b>
<b>Notes</b>	<b>10</b>



# Your Partner in Wound Care

For more information on sizes and ordering codes, please contact  
your local Brennan & Co representative:

<b>Leinster:</b>	Louise French RGN	087 285 4772
<b>Munster:</b>	PJ Mahoney	087 147 8974
<b>Connacht:</b>	Belinda Rushe RGN	087 056 6542
<b>Email:</b>	woundcare@brennanco.ie	01 295 2501



**Brennan & Co**  
Healthcare



# Actico® Cohesive Inelastic Bandage



## Effective for patients

Actico® is a short stretch bandage system, meaning it delivers therapeutic working pressures and tolerable resting pressures.

## Effective for clinicians

*Evidence based:*

Actico® has been tried, tested and loved for over 15 years, generating a wide ranging evidence base, including an RCT.

*Adaptable:*

Actico® offers the solution whether you're treating venous leg ulcers, mixed aetiology ulcers 3,4,5 (ABPI 0.5-0.8)\*, chronic oedema or lymphoedema (with or without leg ulcers). \*with caution – under specialist supervision.

## Safe and easy to apply:

Actico's Safe-Loc® system gives peace of mind for consistently safe, accurate and easy application.

Actico® offers flexible solutions to enable you to treat leg ulcers with or without chronic oedema application.

## Effective for chronic oedema:

Use Actico® with Flexiban® or Cellona® Undercast Padding, with more emphasis on padding for shape improvement during effective reduction in limb size, whilst improving limb proportion. Can be applied as a below knee or full leg bandage.

## Effective for venous leg ulcers:

Use Actico® with Flexiban® or Cellona® Undercast Padding, to protect bony prominences, tendons and fragile tissue, whilst providing effective compression to promote ulcer healing.

## Sizes & Ordering Codes

Size	Order Code
12cm x 6m	88311
10cm x 6m	88310
8cm x 6m	88309
6cm x 6m	88308
4cm x 6m	88307

## Why choose Actico®?

System chosen should...	Consider whether the bandage system...	Actico®
Encourage safe and accurate application.	Is the correct size for the patient's ankle circumference? <b>Is easy to apply?</b> Can be applied accurately and consistently by all those involved in care?	✓
Encourage patient comfort.	Allows the patient to wear appropriate footwear and clothing. Is comfortable and conformable, which encourages patient mobility and confidence.	✓
Be capable of maintaining pressure and be non slip.	Is likely to stay in place, not slip and can be adapted to cope with limb distortion at the next dressing change. Can be adapted to suit limb shape.	✓

# Actico® 2C Two component leg ulcer kit



## Actico® 2C is a two component leg ulcer kit consisting of:

- Layer 1 - a skin friendly comfort layer applied directly to the skin to distribute the compression evenly.
- Layer 2 - an inelastic cohesive compression layer for therapeutic effective compression.
- Latex-free.
- Safe-Loc' system - for consistently safe, accurate and easy application - 'Lock & Roll'.

## Sizes & Ordering Codes

Size	Kit contains	Order Code
Actico®2C - 18-25cm kit	Actico®2C Comfort Layer 1 10cm x 3.5m Actico®2C Compression Layer 2 10cm x 6m	31770
Actico®2C - 25-32cm kit	Actico®2C Comfort Layer 1 10cm x 3.5m + 10cm x 2m Actico®2C Compression Layer 2 10cm x 9m	31771
Actico®2C Comfort Layer 1	10cm x 3.5m	31772
Actico®2C Comfort Layer 2	10cm x 6m	31773

# Actifast® 2 way stretch tubular retention



The versatile ActiFast® 2 way stretch tubular bandage range is easy to apply and can be cut to size and shape.

The ActiFast® range is safe to use where compression is contra-indicated.

## Please note:

When a measurement overlaps, always use the larger size.

The lowest priced 2-way stretch tubular retention bandage on Drug Tariff\*

\*Cost comparison Drug Tariff August 2017.

\*Average reduction of usage in testing against ActiFast® original one way stretch tubular bandage as a result of 'Cut To Fit'.

## Prescriber's guide

Limb/trunk measurement guide:

**Red: Small Limbs**

**Green: Small/Medium Limbs**

**Blue: Large Limbs**

**Yellow: Extra Large Limbs/Heads & Childrens' trunks**

**Beige: Adult Trunks**

**Purple: Extra Large Adult Trunks**

Product	Width (cm)	Length (m)	Order Code
ActiFast® Red	>3.5	10	88574
ActiFast® Green	>5	10	88578
ActiFast® Blue	>7.5	10	88582
ActiFast® Yellow	>10.75	10	88586
ActiFast® Beige	>17.5	10	88588
ActiFast® Purple	>10.75	10	31981

# Activa® Hosiery Liners



Liners provide low compression and are useful for dressing retention and as layer 1 under compression hosiery, as they aid the top layer to glide on more easily.

## Activa® Open Toe/Closed Toe Liners

Colour	Size	Order Code
White Liner - Closed Toe (pack contains x3 liners)		
White	Small	88321
White	Medium	88322
White	Large	88323
White	ExLarge	88324
White	Xtra Xtra Large	88325
Sand Liner - Closed Toe (pack contains x3 liners)		
Sand	Small	88326
Sand	Medium	88327
Sand	Large	88328
Sand	ExLarge	88329
Sand	Xtra Xtra Large	88330
Black Liner - Closed Toe (pack contains x3 liners)		
Black	Small	32840
Black	Medium	32841
Black	Large	32842
Black	ExLarge	32843
Black	Xtra Xtra Large	32844
Sand Liner - <b>Open Toe</b> (pack contains x3 liners)		
Sand	Small	88331
Sand	Medium	88332
Sand	Large	88333
Sand	ExLarge	88334
Sand	Xtra Xtra Large	88335

# Activa® Leg Ulcer Hosiery Kits



Hosiery Kits are an alternative to bandages for healing leg ulcers.

Two garments are worn one on top of the other. A 10mmHg liner followed by higher a compression stocking.

- Optimal 40mmHg\* compression for effective VLU healing.
- Simple to apply consistently and safely.
- Requires minimal training and can be reapplied by the patient, carer and other members of the multi-disciplinary team.

## Activa® Leg Ulcer Hosiery Kit (40mmHg\*)

Colour	Size	Order Code
1 x Sand Stocking BS Class 3 + 2 White Liners		
Sand	Small	88315
Sand	Medium	88316
Sand	Large	88317
Sand	ExLarge	88318
Sand	Xtra Xtra Large	88319
1 x Black Stocking BS Class 3 + 2 Sand Liners		
Black	Small	32675
Black	Medium	32676
Black	Large	32677
Black	ExLarge	32678
Black	Xtra Xtra Large	32679

\*Mean Average

# Flexiban® Padding Bandage



- For protection of exposed bones and nerves under plaster of Paris or synthetic casts.
- Padding underneath splints, zinc paste bandages (Unna's boot), adhesive bandages and compression bandages.
- In wound care: as an airpermeable padding material in combination with textile wound dressings or paraffin gauze dressings.

Size	Order Code	Units Per Pack
10cm x 3.5m	88572	12

# Debrisoft® for debridement



## Debrisoft® Lolly

Ideal for cavity wounds or hard to reach areas (e.g. diabetic foot ulcers, deep pressure ulcers).



## Debrisoft® Pad (10cm x 10cm)

The original Debrisoft® size is suitable for shallow wounds (e.g. venous leg ulcers and other chronic wounds).



## Debrisoft® Pad (13cm x 20cm)

This new bigger size is ideal for larger wounds or areas of hyperkeratosis that may be associated with severe oedema.


## Debrisoft® and Debrisoft® Lolly Ordering Information

Size	Order Code	Units Per Pack
Debrisoft® 5 x 1.9 x 15.5	33224	5
Debrisoft® Pad 10 x 10	34321	5
Debrisoft® Pad 13 x 20	34323	5

# Lomatuell® Pro Gel forming wound contact layer



- Gel formation supports moist wound healing.
- Dressing changes are atraumatic<sup>1</sup>.
- Effective transfer of exudate into the secondary dressing.
- May be used in combination with controlled negative pressure therapy<sup>2</sup>.
- Less disturbance of the wound bed: intervals between dressing changes may be extended.
- Conforms well to the contours of the wound bed.



- A Polymer matrix** elastic, malleable fixing compound, holds petroleum jelly and hydrocolloids in the wide-mesh tulle.
- B Petroleum jelly** moisturises dry wound areas and wound edge.
- C Hydrocolloid** moisturises dry wound areas and wound edge.
- D Open-mesh tulle** carries the polymer matrix and allows transfer of exudate.
- E Gel formation.**

## Sizes & Ordering Codes

**Lomatuell® Pro** Gel forming wound contact net (individually sealed and sterile)

Size (cm)	Order Code	Pack Contains
5 x 5	30870	10
10 x 10	30871	10
10 x 20	30872	10
10 x 30	30873	10

1. LR Data on file, PB FE N-1945/2012

2. Wolber S, et al. International multicentre application study to assess a polyester-tulle primary wound dressing with hydrocolloid particles in terms of usability as well as user and patient satisfaction. Poster presented 24th Conference of the European Wound Management Association (EWMA-GNEAUPP 2014) May 2014 Madrid Spain.

**For more information on sizes and order codes, please contact your local Brennan & Co representative:**

<b>Leinster:</b>	Louise French RGN	087 285 4772
<b>Munster:</b>	PJ Mahoney	087 147 8974
<b>Connacht:</b>	Belinda Rushe RGN	087 056 6542
<b>Email:</b>	woundcare@brennanco.ie	01 295 2501



# ReadyWrap® The easy solution for self-care



**Garments include:** Toe, foot, knee, calf, thigh, extender wrap.

## **ReadyWrap® is better for patients**

Easy to use, with colour coded VELCRO® fasteners.

- Front-fastening short straps aid patient application and readjustment.
- Overlapping straps minimise user application error and can prevent undue loss of pressure if the patient does not apply the product perfectly.
- Low-profile design, allowing patients to wear their usual footwear.
- All products are available in black or beige.

## **ReadyWrap® is clinically effective**

Short stretch properties - Therapeutic working and tolerable lower resting pressure.

- Positive outcomes for leg ulcers, chronic oedema and lymphoedema patients.
- Effective in the management of recurrent swelling (oedema) and ulceration.
- Excellent Static Stiffness Indices, similar to inelastic bandages.
- 50% overlapping straps prevent gapping and can accommodate changing limb circumference to contain and control oedema.

## **ReadyWrap is cost-effective**

ReadyWrap® can provide cost savings associated with:

- Time saved due to quicker and easier application of the device, in comparison to compression bandages.
- Reduced use of bandaging materials – ReadyWrap can be washed and reused over a period of six months.
- Reducing clinic visits due to self-care at home.
- Educating nursing time as it can be easily applied by the patient or their carer.

**For more information on sizes and order codes, please contact your local Brennan & Co representative:**

<b>Leinster:</b>	Louise French RGN	087 285 4772
<b>Munster:</b>	PJ Mahoney	087 147 8974
<b>Connacht:</b>	Belinda Rushe RGN	087 056 6542
<b>Email:</b>	woundcare@brennanco.ie	01 295 2501

# Vliwasorb Adhesive® superabsorbent wound dressing compression



Adhesive superabsorbent wound dressing with moderate to high levels of exudate.


## Key benefits:

- Effective exudate management.
- Convenience and comfort of a skin-friendly, self-adjusting adhesive.

## Ordering Information

Vliwasorb Adhesive Superabsorbent (10 singles)  
Individually sealed and sterile

Size	Order Code
12cm x 12cm	30991
15cm x 15cm	30992
15cm x 25cm	30993



**Flexible wound contact layer** prevents adhesion to the wound and means a primary contact layer is not required.

**Distribution layer** exudate is quickly and evenly distributed into the absorbent core.

**Absorbent core with superabsorbent polymers (SAP)** binds large amounts of wound exudate - protecting against maceration.

**Blue nonwoven clothing protection** protects against soiling.

# Vliwasorb® Pro Highly conformable dressing for highly exuding wounds



## Did you know?

Vliwasorb® Pro may be used as a primary wound dressing - it's integrated wound contact layer negates the need to use an additional contact layer underneath - saving budget.

## The structure of Vliwasorb® Pro

Non-wound side

Wound side

Rounded edges and ergonomic shape **for good conformability.**

Wide edges **for high wearing comfort.**

Blue clothing protection layer **protects clothes and bed linen.**

Superabsorbent core with distribution layer **for effective exudate management and flexibility.**

Smooth wound contact layer **for atraumatic dressing changes.**

## Ordering Information

Individually sealed and sterile

Size (cm)	Wound pad (cm)	Order Code
12.5 x 12.5	8.5 x 8.5	32641
12.5 x 22.5	8.5 x 18.5	32642
22 x 22	18 x 18	32643
22 x 32	18 x 28	32644

## Notes

[illegible]



**For more information on sizes and order codes, please contact your local Brennan & Co representative:**

<b>Leinster:</b>	Louise French RGN	087 285 4772
<b>Munster:</b>	PJ Mahoney	087 147 8974
<b>Connacht:</b>	Belinda Rushe RGN	087 056 6542
<b>Email:</b>	woundcare@brennanco.ie	01 295 2501